



www.pro-lab.com

CERTIFICATE OF ANALYSIS

Product Name: SPOT INDOLE REAGENT	Lot Number: K32274
Product Code: PL.391-10	Expiration Date: 2022 09 30

Additional Information

Description: Pale orange solution.	Concentrate Lot number: n/a Production date: 2021/09/06	Storage conditions: 2-30°C Precautions: refer to SDS
------------------------------------	--	---


Quality Control Testing

Challenged with test organisms:
Escherichia coli NCTC 12241
Neisseria gonorrhoea NCTC 12700

Test results:
 Positive reaction: the development of a blue colour within 3 minutes with *E.coli*.
 Negative reaction: the development of a pink colour within 3 minutes with *N.gonorrhoea* .

Comments:
 We hereby certify that samples of the above Lot have been tested in accordance with current specifications and that the above is a true record of the results obtained.

Release

2021/09/15	
Date	Quality Control

Canada ☎ 800 268 2341 Fax 905 731 0206
 20 Mural Street, Unit #4, Richmond Hill, Ontario,
 L4B 1K3
 support@pro-lab.com

U.S.A. ☎ (512) 832 9145 Fax (512) 832 6424
 21 Cypress Blvd., Suite 1155, Round Rock,
 Texas 78665-1046
 support@pro-lab.us

U.K. ☎ 0151 353 1613 Fax 0151 353 1614
 3 Bassendale Road, Bromborough, Wirral,
 Merseyside CH62 3QL
 uksupport@pro-lab.com