

### CERTIFICATE OF ANALYSIS

<b>Product Name:</b> Mucolyse	<b>Lot Number:</b> C12946	
<b>Product Code:</b> PL.701	<b>Expiration Date:</b> 2017 07	

#### Additional Information


Storage Conditions:	2-8°C
Precautions:	Refer to Instructions for Use

#### Quality Control Testing

	Criteria	Result
<b>Appearance</b>		
Pellet	Compact white pellet	PASS
Reconstituted Supplement	Clear colourless solution	PASS
<b>Sterility Test</b>		
	Conforms	PASS
<b>Organism</b>		
<i>Haemophilus influenzae</i> ATCC 10211	<b>Growth at 24 hours*</b> Good Growth	PASS
<i>Streptococcus pneumoniae</i> ATCC 6305	Good Growth	PASS

\*Compared to control plates not containing supplement.

#### Release

2015 08 17	
Date	Quality Control