



www.pro-lab.com

CERTIFICATE OF ANALYSIS

Product Name: TB Staining Kit (Malachite Green)	Lot Number: K17466
Product Code: PL.8061/25	Expiration Date: 2014 12

Product Code	Quantity	Product description
PL.7018/25	1 X 250ml	ZN Carbol Fuchsin
PL.7024/25	2 X 250ml	Differentiator for ZN & Kinyoun Carbol Fuchsin
PL.7030/25	1 X 250ml	Malachite Green

Additional Information

Description: Magenta, clear and green solutions.	Concentrate Lot numbers: 5896/5933 Production date: 2013/12/23	Storage conditions: 15-25°C Precautions: refer to MSDS
--	---	---

Quality Control Testing

Challenged with test organisms:
Mycobacterium tuberculosis H37 Rv NCTC 7416
Streptomyces griseus NCTC 7807
Cryptosporidia var parvum wild strain

Test results:
 Conforms to performance specifications when used in standard Ziehl Neelson stain method, giving crisp red staining of Acid/Alcohol Fast bacilli with a clean counterstained background.

Comments:
 We hereby certify that samples of the above Lot have been tested in accordance with current specifications and that the above is a true record of the results obtained.

Release

13/12/23	<i>maun .</i>
Date	Quality Control

Canada ☎ 800 268 2341 Fax 905 731 0206
 20 Mural Street, Unit #4, Richmond Hill, Ontario,
 support@pro-lab.com

U.S.A. ☎ 800 522 7740 Fax 800 332 0450
 21 Cypress Blvd., Suite 1070, Round Rock,
 ussupport@pro-lab.com

U.K. ☎ 0151 353 1613 Fax 0151 353 1614
 3 Bassendale Road, Bromborough, Wirral,
 uksupport@pro-lab.com