




www.pro-lab.com

CERTIFICATE OF ANALYSIS

Product Name: TB Staining Kit (Malachite Green)		Lot Number: K31687
Product Code: PL.8061/25		Expiration Date: 2022 05 31
Product Code	Quantity	Product description
PL.7018/25	1 X 250ml	ZN Carbol Fuchsin
PL.7024/25	2 X 250ml	Differentiator for ZN & Kinyoun Carbol Fuchsin
PL.7030/25	1 X 250ml	Malachite Green
Additional Information		
Description: Magenta, clear and green solutions.	Concentrate Lot numbers: K31628/K31641 Production date: 2021/05/17	Storage conditions: 15-25°C Precautions: refer to SDS
Quality Control Testing		
<p><u>Challenged with test organisms:</u> <i>Mycobacterium scrofulaceum</i> NCTC 10803 <i>Escherichia coli</i> NCTC 12241</p> <p><u>Test results:</u> Conforms to performance specifications when used in standard Ziehl Neelson stain method, giving crisp red staining of Acid/Alcohol Fast bacilli with a clean counterstained background.</p> <p><u>Comments:</u> We hereby certify that samples of the above Lot have been tested in accordance with current specifications and that the above is a true record of the results obtained.</p>		
Release		
2021/05/21		
Date	Quality Control	

Canada ☎ 800 268 2341 Fax 905 731 0206
 20 Mural Street, Unit #4, Richmond Hill,
 Ontario, L4B 1K3
 support@pro-lab.com

U.S.A. ☎ (512) 832 9145 Fax (512) 832 6424
 21 Cypress Blvd., Suite 1155, Round Rock,
 Texas 78665-1046
 support@pro-lab.us

U.K. ☎ 0151 353 1613 Fax 0151 353 1614
 3 Bassendale Road, Bromborough,
 Wirral, Merseyside, CH62 3QL
 uksupport@pro-lab.com