

## QUALITY CONTROL CERTIFICATE

### **HIPPURATE HYDROLYSIS Diagnostic Tablets**

REF: 56711  
LOT: 1301-1  
Date of analysis: 06-2013  
Date of expiry: 06-2017  
Appearance: Plane tablet, 6 mm in diameter.  
Complies.  
Average weight: 90 mg.  
Complies.

#### Performance:

Test strains:

*Streptococcus agalactiae* ATCC 12386  
*Streptococcus pyogenes* ATCC 12344  
*Staphylococcus aureus* 25923  
*E. coli* ATCC 25922

Product complies with specification.

ROSCO DIAGNOSTICA A/S



José B. Casals  
Quality Assurance Manager