

## QUALITY CONTROL CERTIFICATE

### **VOGES-PROSKAUER Diagnostic Tablets**

REF: 57711  
LOT: 1701-1  
Date of analysis: 2017-08  
Date of expiry: 2020-02  
Appearance: Plane tablet, 6 mm in diameter.  
Complies.  
Average weight: 90 mg.  
Complies.


#### Performance:

Test strains:

*E. coli* ATCC 25922  
*Staphylococcus aureus* ATCC 25923  
*Enterobacter cloacae* ATCC 13047

Product complies with specification.

Rosco Diagnostica A/S



Michela Venturi  
Manager Quality Assurance