

## QUALITY CONTROL CERTIFICATE

### **Rapid CARB Sreen Kit Diagnostic Tablets**

REF: 98021  
LOT: 1301-2  
Date of analysis: 05-2013, 06-2013  
Date of expiry: 05-2016  
Appearance: Plane tablet, 6 mm in diameter.  
Complies.  
Average weight: 90 mg.  
Complies.

### Performance:

#### Test strains

DIATABS	Imipenem + Indicator (CARB)	Negative Control (CARB)
<i>Klebsiella pneumoniae</i> BAA 1705	Positive	Negative
<i>E. coli</i> ATCC 25922	Negative	Negative

Product complies with specifications.

ROSCO DIAGNOSTICA A/S



José B. Casals  
Quality Assurance Manager