

Certificate of Quality Control

Manufacturer **SIFIN** Institut für Immunpräparate und Nährmedien GmbH Berlin,
Berliner Allee 317-321, 13088 Berlin, Germany

Product **Enteroclon
Anti-Salmonella I (A – E)**

Catalogue No. (Package size) TR 1111 (1 ml)

Lot No. **116 08 12**

Expiry date **2015-08**

Intended use Enteroclon Anti-Salmonella I (A – E) records all Salmonella serovars of the groups A to E in slide agglutination.

Specificity No agglutination with Salmonella of the groups F to 67 as well as Shigella sonnei (S and F form) and all serotypes of Shigella flexneri.

Preservatives Sodium azide 0.9 mg/ ml

Storage 2...8 °C

Berlin,

27.09.13


Head of Quality Control