

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
 ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh


Berliner Allee 317-321
 13088 Berlin, Germany

T: +49 30 927 030-0
 F: +49 30 927 030-30
 info@sifin.de
 www.sifin.de

Product	Anti-Salmonella I (A-E)
Catalogue no.	TR 1111
Package size	1 ml
Lot no.	154 08 14
Expiry date	2016-08
Intended use	Anti-Salmonella I (A-E) records all <i>Salmonella</i> serovars of the groups A-E and the Vi-antigen in slide agglutination.
Specificity	No agglutination with <i>Salmonella</i> strains of the group F-67 as well as <i>Shigella sonnei</i> and all serotypes of <i>Shigella flexneri</i> .
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin, 16.10.14

 place, date



 signature Head of Quality Control