



Certificate of Quality Control

Manufacturer

SIFIN Institut für Immunpräparate und Nährmedien GmbH Berlin,

Berliner Allee 317-321, 13088 Berlin, Germany

Product

Enterocion

Anti-Salmonella II (F - 67)

Catalogue No. (Package size)

TR 1121 (1 ml)

Lot No.

170 11 12

Expiry date

2016-11

Intended use

Enteroclon Anti-Salmonella II (F-67) records all Salmonella serovars of the groups F to 67 in slide agglutination.

Specificity

No agglutination with Salmonella of the groups A to E as well as Shigella sonnei (S and F form) and all serotypes of Shigella flexneri. Due to related antigens among genera of the Enterobacteriaceae family there is the possibility of agglutination reactions, for instance, with strains of Proteus vulgaris and Proteus mirabilis as well as representatives of Citrobacter species.

Preservatives

Sodium azide 0.9 mg/ ml

Storage

2...8 °C

Berlin,