

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321 13088 Berlin, Germany

T: +49 30 927 030-0 F: +49 30 927 030-30 info@sifin.de www.sifin.de

Product

Anti-Salmonella II (F-67)

Catalogue no.

TR 1121

Package size

1 ml

Lot no.

7860917

Expiry date

2021-09

Intended use

Anti-Salmonella II (F-67) records all Salmonella strains of the groups F-67

in slide agglutination.

Specificity

No agglutination with Salmonella of the groups A to E as well as Shigella sonnei. Due to related antigens among genera of the Enterobacteriaceae family there is the possibility of agglutination reactions, for instance, with strains of Proteus vulgaris and Proteus mirabilis as well as representatives

of Citrobacter.

Preservatives

Sodium azide 0.9 mg/ml

Storage

2 °C...8 °C

Berlin,

place, date

13.10.2217

signature Head of Quality Control

T. Solvey