

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321
13088 Berlin, Germany

T: +49 30 700 144-0
F: +49 30 700 144-30
info@sifin.de
www.sifin.de

Product	Anti-Salmonella II (F-67)
Catalogue no.	TR1121
Package size	1 ml
Lot no.	4340621
Expiry date	2025-05
Intended use	Anti-Salmonella II (F-67) records all <i>Salmonella</i> strains of the groups F-67 in slide agglutination.
Specificity	No agglutination with <i>Salmonella</i> of the groups A to E as well as <i>Shigella sonnei</i> . Due to related antigens among genera of the <i>Enterobacteriaceae</i> family there is the possibility of agglutination reactions, for instance, with strains of <i>Proteus vulgaris</i> and <i>Proteus mirabilis</i> as well as representatives of <i>Citrobacter</i> .
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin,

23.06.21

place, date

J. Lueders

signature Head of Quality Control