

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
 ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321
 13088 Berlin, Germany

T: +49 30 927 030-0
 F: +49 30 927 030-30
 info@sifin.de
 www.sifin.de

Product	Anti-Salmonella OMB O-Group-Pool: C, F, G, H
Catalogue no.	TR 1161
Package size	1 ml
9091018	9091018
Expiry date	2020-10
Intended use	Anti-Salmonella OMB records all <i>Salmonella</i> strains belonging to the O-Groups C, F, G, H in slide agglutination.
Specificity	No agglutination with <i>Salmonella</i> strains which are recorded by Anti-Salmonella OMA (O-Group-Pool: A, B, D, E, L), Anti-Salmonella OMC (O-Group-Pool: I, J, K, M, N, O, P), Anti-Salmonella OMD (O-Group-Pool: Q, R, S, T, U, V, W), Anti-Salmonella OME (O-Group-Pool: X, Y, Z, 51, 52, 53), Anti-Salmonella OMF (O-Group-Pool: 54-59) and Anti-Salmonella OMG (O-Group-Pool: 60-63, 65-67).
Preservatives	Sodium azide 0,9 mg/ml
Storage	2 °C... 8 °C

Berlin,

 place, date

8.11.2018



 signature Head of Quality Control