

Qualitätskontrollzertifikat

Manufacturer	SIFIN Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin
Product	Test Reagent Anti-Salmonella OMG (O-Group-Pool: 60-63, 65-67)
Catalogue No. (Package size)	TR 1174 (1 ml)
Lot No.	64 05 13
Expiry date	2015-05
Intended use	Test Reagent Anti-Salmonella OMG records all <i>Salmonella</i> strains belonging to the O-Groups 60-63, 65-67.
Specificity	No agglutination with <i>Salmonella</i> strains which are recorded by Test Reagents Anti-Salmonella OMA (O-Group-Pool: A, B, D, E, L), Anti-Salmonella OMB (O-Group-Pool: C, F, G, H), Anti-Salmonella OMC (O-Group-Pool: I, J, K, M, N, O, P), Anti-Salmonella OMD (O-Group-Pool: Q, R, S, T, U, V, W), Anti-Salmonella OME (O-Group-Pool: X, Y, Z, 51, 52, 53) and Anti- Salmonella OMF (O-Group-Pool: 54-59).
Preservatives	Sodium azide 0,9 mg/ml
Storage	2...8 °C

Berlin, 31.05.13

J. Kerschling
Head of Quality Control