

## Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321  
13088 Berlin, Germany

T: +49 30 927 030-0  
F: +49 30 927 030-30  
info@sifin.de  
www.sifin.de

<b>Product</b>	<b>Anti-Salmonella Group C</b>
Catalogue no.	TR 1202
Package size	1 ml
<b>Lot no.</b>	<b>0861217</b>
<b>Expiry date</b>	<b>2020-12</b>
Intended use	Anti-Salmonella Group C records all <i>Salmonella</i> strains of group C (C <sub>1</sub> ; C <sub>2</sub> -C <sub>3</sub> ) in slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> O antigens and <i>Salmonella</i> H antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2°C...8 °C

Berlin,  
place, date

20.01.2018

  
signature Head of Quality Control