

Certificate of Quality Control

Manufacturer **SIFIN** Institut für Immunpräparate und Nährmedien GmbH Berlin,
Berliner Allee 317-321, 13088 Berlin, Germany

Product **Enteroclon Anti-Salmonella O 7**

Catalogue No. TR 1305 (1 ml)
(Package size)

Lot No. **83 05 12**

Expiry date **2015-05**

Intended use Enteroclon Anti-Salmonella O 7 records all Salmonella strains as
long as they contain antigen O 7.

Specificity No agglutination with heterologous Salmonella O antigens and
Salmonella H antigens.

Preservatives Sodium azide 0.9 mg / ml

Storage 2...8 °C

Berlin, 11.09.12


Head of Quality Control