


Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

Product	Anti-Salmonella O:7
Catalogue no.	TR1305
Package size	1 ml
Lot no.	3380421
Expiry date	2024-03
Intended use	Anti-Salmonella O:7 records all <i>Salmonella</i> strains, as long as they contain antigen O:7 in slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> O antigens and <i>Salmonella</i> H antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C... 8 °C

Berlin, 21.06.21
place, date


signature Head of Quality Control