


Certificate of Quality Control

Manufacturer	SIFIN Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany
Product	Enteroclon Anti-Salmonella O 8
Catalogue No. (Package size)	TR 1306 (1 ml)
Lot No.	112 10 09
Expiry date	2014-11
Intended use	Enteroclon Anti-Salmonella O 8 records all Salmonella strains as long as they contain antigen O 8 (group C ₂ -C ₃).
Specificity	No agglutination with heterologous Salmonella O antigens and Salmonella H antigens.
Preservatives	Sodium azide 0.9 mg / ml
Storage	2...8 °C

Berlin, 21.12.08



Head of Quality Control