

Certificate of Quality Control

Manufacturer **SIFIN** Institut für Immunpräparate und Nährmedien GmbH Berlin,
Berliner Allee 317-321, 13088 Berlin, Germany

Product **Enteroclon Anti-Salmonella O 8**

Catalogue No.
(Package size) TR 1306 (1 ml)

Lot No. **189 11 13**

Expiry date **2018-11**

Intended use Enteroclon Anti-Salmonella O 8 records all Salmonella strains as
long as they contain antigen O 8 (group C₂-C₃).

Specificity No agglutination with heterologous Salmonella O antigens and
Salmonella H antigens.

Preservatives Sodium azide 0.9 mg / ml

Storage 2...8 °C

Berlin,

15.07.14


Head of Quality Control