

Certificate of Quality Control

Manufacturer	SIFIN Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany
Product	Enteroclon Anti-Salmonella O 14
Catalogue No. (Package size)	TR 1309 (1 ml)
Lot No.	51 05 13
Expiry date	2015-05
Intended use	Enteroclon Anti-Salmonella O 14 records all Salmonella strains, as long as they contain of group H (6,14).
Specificity	No agglutination with heterologous Salmonella O antigens and Salmonella H antigens
Preservatives	Sodium azide 0.9 mg / ml
Storage	2...8 °C

Berlin,

23.05.73

F. Ludwig

Head of Quality Control