


Certificate of Quality Control

Manufacturer	SIFIN Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany
Product	Enteroclon Anti-Salmonella O 15
Catalogue No. (Package size)	TR 1310 (1 ml)
Lot No.	23 02 12
Expiry date	2016-02
Intended use	Enteroclon Anti-Salmonella O 15 records all Salmonella strains of group E1 after conversion by the phages ϵ_{15} or ϵ_{15} and ϵ_{34} . Antigens O 15 or O 15,34 replace antigen O 10. Antigen O 10 is no longer detectable after lysogenisation.
Specificity	No other agglutinations with heterologous Salmonella O antigens and Salmonella H antigens are known.
Preservatives	Sodium azide 0.9 mg / ml
Storage	2...8 °C

Berlin,

17.04.12



Head of Quality Control