

## Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321  
13088 Berlin, Germany

T: +49 30 700 144-0  
F: +49 30 700 144-30  
info@sifin.de  
www.sifin.de

<b>Product</b>	<b>Anti-Salmonella O:15</b>
Catalogue no.	TR1310
Package size	1 ml
<b>Lot no.</b>	<b>7410920</b>
<b>Expiry date</b>	<b>2024-08</b>
Intended use	Anti-Salmonella O:15 records all <i>Salmonella</i> strains of group E1 after conversion by the phages $\epsilon$ 15 or $\epsilon$ 15 and $\epsilon$ 34. Antigens O:15 or O:15,34 replace antigen O:10. Antigen O:10 is no longer detectable after lysogenisation.
Specificity	No agglutination with heterologous <i>Salmonella</i> O antigens and <i>Salmonella</i> H antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin,

07.10.2020

place, date

*F. Lauer*

signature Head of Quality Control