

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321
13088 Berlin, Germany

T: +49 30 927 030-0
F: +49 30 927 030-30
info@sifin.de
www.sifin.de

Product	Anti-Salmonella O:19
Catalogue no.	TR 1311
Package size	1 ml
Lot no.	9571117
Expiry date	2022-11
Intended use	Anti-Salmonella O:19 records all <i>Salmonella</i> strains of group E ₄ (O:1,3,19).
Specificity	Anti-Salmonella O:19 cross-reacts slightly with several strains of group B with O:27. However, since the test reagent is used only once the <i>Salmonella</i> strain's affiliation to group E has been confirmed, this cross-reaction is irrelevant to the serotyping procedure. No other agglutinations with heterologous <i>Salmonella</i> O antigens and <i>Salmonella</i> H antigens are known.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin, 21.11.2017
place, date

F. Herdewig
signature Head of Quality Control