



## Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485


sifin diagnostics gmbh

Berliner Allee 317-321  
13088 Berlin, Germany

T: +49 30 927 030-0  
F: +49 30 927 030-30  
info@sifin.de  
www.sifin.de

<b>Product</b>	<b>Anti-Salmonella O:20</b>
Catalogue no.	TR 1312
Package size	1 ml
<b>Lot no.</b>	<b>16 01 15</b>
<b>Expiry date</b>	<b>2019-01</b>
Intended use	Anti-Salmonella O:20 records all <i>Salmonella</i> strains of group C <sub>2</sub> -C <sub>3</sub> , as long as they contain antigen O:20 in slide agglutination
Specificity	No agglutination with heterologous <i>Salmonella</i> O antigens and <i>Salmonella</i> H antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin, 18.02.15  
\_\_\_\_\_  
place, date

  
\_\_\_\_\_  
signature Head of Quality Control