

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

Product	Anti-Salmonella O:20
Catalogue no.	TR 1312
Package size	1 ml
Lot no.	7240819
Expiry date	2023-08
Intended use	Anti-Salmonella O:20 records all <i>Salmonella</i> strains of group C ₂ -C ₃ , as long as they contain antigen O:20 in slide agglutination
Specificity	No agglutination with heterologous <i>Salmonella</i> O antigens and <i>Salmonella</i> H antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C... 8 °C

Berlin,

13.09.2018

place, date



signature Head of Quality Control