

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

Product	Anti-Salmonella O:34
Catalogue no.	TR 1314
Package size	1 ml
Lot no.	0131119
Expiry date	2024-11
Intended use	Anti-Salmonella O:34 records all <i>Salmonella</i> strains of group E ₁ , as long as they contain antigen O:34 in slide agglutination.
Specificity	Anti-Salmonella O:34 shows cross reactions with O:12 ₂ , O:28 and O:55 and should therefore used only for the differentiation in group E. No agglutinations with <i>Salmonella</i> H antigens are known.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin,

28.01.2020

place, date

F. Leuch

signature Head of Quality Control