

Certificate of Quality Control

Manufacturer **SIFIN** Institut für Immunpräparate und Nährmedien GmbH Berlin,
Berliner Allee 317-321, 13088 Berlin

Product **Enteroclon Anti-Salmonella O 46**

Catalogue No.
(Package size) TR 1315 (1 ml)

Lot No. **79 06 13**

Expiry date **2017-06**

Intended use Enteroclon Anti-Salmonella O 46 records all Salmonella strains of group D₂ (O 9,46) in the slide agglutination test. Strains of group D₃ (O 9,12,46,27) do not agglutinate or agglutinate only weakly.

Specificity No agglutination with heterologous Salmonella O antigens and Salmonella H antigens.

Preservatives Sodium azide 0.9 mg/ ml

Storage 2...8 °C

Berlin,

15.08.13


Head of Quality Control