

Certificate of Quality Control

Manufacturer	SIFIN Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany
Product	Enteroclon Anti-Salmonella O 39
Catalogue No. (Package size)	TR 1345 (1 ml)
Lot No.	142 10 11
Expiry date	2016-10
Intended use	Enteroclon Anti-Salmonella O 39 records all Salmonella strains, as long as they contain antigen O 39.
Specificity	No other agglutinations with heterologous Salmonella O antigens and Salmonella H antigens are known.
Preservatives	Sodium azide 0.9 mg / ml
Storage	2...8 °C

Berlin, 21.11.11

K. P. G. J.
Head of Quality Control