


---

## Certificate of Quality Control

Manufacturer	<b>SIFIN</b> Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany
Product	<b>Enteroclon Anti-Salmonella O 43</b>
Catalogue No. (Package size)	TR 1349 (1 ml)
Lot No.	<b>148 11 11</b>
Expiry date	<b>2016-11</b>
Intended use	Enteroclon Anti-Salmonella O 43 records all Salmonella strains, as long as they contain antigen O 43 (Group O 43).
Specificity	No other agglutinations with heterologous Salmonella O antigens and Salmonella H antigens are known.
Preservatives	Sodium azide 0.9 mg / ml
Storage	2...8 °C

Berlin, 22.12.11

  
\_\_\_\_\_  
Head of Quality Control