

## Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321  
13088 Berlin, Germany

T: +49 30 927 030-0  
F: +49 30 927 030-30  
info@sifin.de  
www.sifin.de

<b>Product</b>	<b>Anti-Salmonella O:47</b>
Catalogue no.	TR 1353
Package size	1 ml
<b>Lot no.</b>	<b>44 03 13L</b>
<b>Expiry date</b>	<b>2018-03</b>
Intended use	Anti-Salmonella O:47 records all <i>Salmonella</i> strains, as long as they contain antigen O:47 in slide agglutination
Specificity	No agglutination with heterologous <i>Salmonella</i> O antigens and <i>Salmonella</i> H antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin, 04.04.14  
\_\_\_\_\_  
place, date

  
\_\_\_\_\_  
signature Head of Quality Control