

Certificate of Quality Control

Manufacturer	SIFIN Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany
Product	Enteroclon Anti-Salmonella O 52
Catalogue No. (Package size)	TR 1357 (1 ml)
Lot No.	30 02 13
Expiry date	2018-02
Intended use	Enteroclon Anti-Salmonella O 52 records all Salmonella strains, as long as they contain antigen O 52.
Specificity	No agglutination with heterologous Salmonella O antigens and Salmonella H antigens.
Preservatives	Sodium azide 0.9 mg/ ml
Storage	2...8 °C

Berlin,

05.08.13



Head of Quality Control