

### Certificate of Quality Control

Manufacturer	<b>SIFIN</b> Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany
Product	<b>Enteroclon Anti-Salmonella O 55</b>
Catalogue No. (Package size)	TR 1360 (1 ml)
Lot No.	<b>105 06 12</b>
Expiry date	<b>2017-06</b>
Intended use	Enteroclon Anti-Salmonella O 55 records all Salmonella strains, as long as they contain antigen O 55.
Specificity	No agglutination with heterologous Salmonella O antigens and Salmonella H antigens, with the exception of antigen O 59 with some strains of this group.
Preservatives	Sodium azide 0.9 mg/ ml
Storage	2...8 °C

Berlin, 27.07.12

  
Head of Quality Control