


---

### Certificate of Quality Control

Manufacturer	<b>SIFIN</b> Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany
Product	<b>Enteroclon Anti-Salmonella O 65</b>
Catalogue No. (Package size)	TR 1369 (1 ml)
Lot No.	<b>24 02 12</b>
Expiry date	<b>2017-02</b>
Intended use	Enteroclon Anti-Salmonella O 65 records all Salmonella strains, as long as they contain antigen O 65.
Specificity	No agglutination with heterologous Salmonella O antigens and Salmonella H antigens.
Preservatives	Sodium azide 0.9 mg/ ml
Storage	2...8 °C

Berlin,

17.04.12

  
Head of Quality Control