

Certificate of Quality Control

Manufacturer	SIFIN Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany
Product	Enteroclon Anti-Salmonella H b
Catalogue No. (Package size)	TR 1402 (1 ml)
Lot No.	143 09 12
Expiry date	2015-09
Intended use	Enteroclon Anti-Salmonella H b records all Salmonella strains, as long as they contain antigen H b.
Specificity	No agglutination with heterologous Salmonella H antigens and Salmonella O antigens
Preservatives	Sodium azide 0.9 mg / ml
Storage	2...8 °C

Berlin, 19.12.12

K. Reyer
Head of Quality Control