

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
 ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321
 13088 Berlin, Germany

T: +49 30 927 030-0
 F: +49 30 927 030-30
 info@sifin.de
 www.sifin.de

| | |
|--------------------|--|
| Product | Anti-Salmonella H:b |
| Catalogue no. | TR 1402 |
| Package size | 1 ml |
| Lot no. | 7620817 |
| Expiry date | 2022-08 |
| Intended use | Anti-Salmonella H:b records all <i>Salmonella</i> strains, as long as they contain antigen H:b in slide agglutination. |
| Specificity | No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens. |
| Preservatives | Sodium azide 0.9 mg/ml |
| Storage | 2 °C...8 °C |

Berlin,

20.09.2017

place, date



signature Head of Quality Control