

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

Product	Anti-Salmonella H:b
Catalogue no.	TR1402
Package size	1 ml
Lot no.	6350820
Expiry date	2023-07
Intended use	Anti-Salmonella H:b records all <i>Salmonella</i> strains, as long as they contain antigen H:b in slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin,

24. 09. 2020

place, date

F. Ludwig

signature Head of Quality Control