

## Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
 ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321  
 13088 Berlin, Germany

T: +49 30 927 030-0  
 F: +49 30 927 030-30  
 info@sifin.de  
 www.sifin.de

<b>Product</b>	<b>Anti-Salmonella H:d</b>
Catalogue no.	TR 1404
Package size	1 ml
<b>Lot no.</b>	<b>46 03 15</b>
<b>Expiry date</b>	<b>2018-04</b>
Intended use	Anti-Salmonella H:d records all <i>Salmonella</i> strains, as long as they contain antigen H:d in slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C... 8 °C

Berlin, 23.04.15  
 place, date

  
 signature Head of Quality Control