

Certificate of Quality Control

Manufacturer	SIFIN Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany
Product	Enteroclon Anti-Salmonella H d
Catalogue No. (Package size)	TR 1404 (1 ml)
Lot No.	93 06 13
Expiry date	2016-06
Intended use	Enteroclon Anti-Salmonella H d records all Salmonella strains, as long as they contain antigen H d in slide agglutination.
Specificity	No agglutination with heterologous Salmonella H antigens and Salmonella O antigens.
Preservatives	Sodium azide 0.9 mg/ ml
Storage	2...8 °C

Berlin,

05.08.13


Head of Quality Control