

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
 ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321
 13088 Berlin, Germany

T: +49 30 927 030-0
 F: +49 30 927 030-30
 info@sifin.de
 www.sifin.de

Product	Anti-Salmonella H:E
Catalogue no.	TR 1405
Package size	1 ml
Lot no.	65 04 16
Expiry date	2019-04
Intended use	Anti-Salmonella H:E records all <i>Salmonella</i> strains, as long as they contain antigen H:e (e.g. H:e,h; H:e,n,x; H:e,n,z ₁₅ ; H:e,n,x,z ₁₅) in slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin, 27.05.2016
 place, date


 signature Head of Quality Control