

## Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

<b>Product</b>	<b>Anti-Salmonella H:E</b>
Catalogue no.	TR1405
Package size	1 ml
<b>Lot no.</b>	<b>2530320</b>
<b>Expiry date</b>	<b>2023-02</b>
Intended use	Anti-Salmonella H:E records all <i>Salmonella</i> strains, as long as they contain antigen H:e (e.g. H:e,h; H:e,n,x; H:e,n,z <sub>15</sub> ; H:e,n,x,z <sub>15</sub> ) in slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin, 01.09.2020  
\_\_\_\_\_  
place, date

7. Ludwig  
\_\_\_\_\_  
signature Head of Quality Control