

# sifin

## Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485


sifin diagnostics gmbh

Berliner Allee 317-321  
13088 Berlin, Germany

T: +49 30 927 030-0  
F: +49 30 927 030-30  
info@sifin.de  
www.sifin.de

<b>Product</b>	<b>Anti-Salmonella H:f</b>
Catalogue no.	TR 1407
Package size	1 ml
<b>Lot no.</b>	<b>132 07 14</b>
<b>Expiry date</b>	<b>2017-07</b>
Intended use	Anti-Salmonella H:f records all <i>Salmonella</i> strains, as long as they contain in an antigen combination the H f antigen (e. g. H: f,g; H:f,g,s; H:f,g,t; H:f,g,m,t). In slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin, 09 09 14  
\_\_\_\_\_  
place, date

  
\_\_\_\_\_  
signature Head of Quality Control