

## Certificate of Quality Control

Manufacturer	<b>SIFIN</b> Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany
Product	<b>Enteroclon Anti-Salmonella H f</b>
Catalogue No. (Package size)	TR 1407 (1 ml)
Lot No.	<b>189 12 12</b>
Expiry date	<b>2015-12</b>
Intended use	Enteroclon Anti-Salmonella H f records all Salmonella strains, as long as they contain in an antigen combination the H f antigen (e. g. f,g; f,g,s; f,g,t; f,g,m,t).
Specificity	No agglutination with heterologous Salmonella H antigens and Salmonella O antigens
Preservatives	Sodium azide 0.9 mg / ml
Storage	2...8 °C

Berlin, 10.01.13

K. Pflüger  
Head of Quality Control