

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
 ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321
 13088 Berlin, Germany

T: +49 30 927 030-0
 F: +49 30 927 030-30
 info@sifin.de
 www.sifin.de

Product	Anti-Salmonella H:f
Catalogue no.	TR 1407
Package size	1 ml
Lot no.	1030119
Expiry date	2022-01
Intended use	Anti-Salmonella H:f records all <i>Salmonella</i> strains, as long as they contain in an antigen combination the H f antigen (e. g. H:f,g; H:f,g,s; H:f,g,t; H:f,g,m,t). In slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin, 14.03.2019
 place, date

7. Burdwig
 signature Head of Quality Control