

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321
13088 Berlin, Germany

T: +49 30 700 144-0
F: +49 30 700 144-30
info@sifin.de
www.sifin.de

Product	Anti-Salmonella H:f
Catalogue no.	TR1407
Package size	1 ml
Lot no.	1830220
Expiry date	2023-01
Intended use	Anti-Salmonella H:f records all <i>Salmonella</i> strains, as long as they contain in an antigen combination the H f antigen (e. g. H:f,g; H:f,g,s; H:f,g,t; H:f,g,m,t). In slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin,

place, date

08.09.2020



signature Head of Quality Control