

Certificate of Quality Control

Manufacturer	SIFIN Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany
Product	Enteroclon Anti-Salmonella H h
Catalogue No. (Package size)	TR 1409 (1 ml)
Lot No.	101 06 12
Expiry date	2015-06
Intended use	Enteroclon Anti-Salmonella H h records all Salmonella strains, as long as they contain antigen H h (in the antigen combination e,h) in slide agglutination.
Specificity	no agglutination with Salmonella O antigens and heterologous Salmonella H antigens
Preservatives	Sodium azide 0.9 mg/ ml
Storage	2...8 °C

Berlin, 20.07.12


Head of Quality Control