

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
 ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh


Berliner Allee 317-321
 13088 Berlin, Germany

T: +49 30 927 030-0
 F: +49 30 927 030-30
 info@sifin.de
 www.sifin.de

Product	Anti-Salmonella H:h
Catalogue no.	TR 1409
Package size	1 ml
Lot no.	250 12 15
Expiry date	2018-12
Intended use	Anti-Salmonella H:h records all <i>Salmonella</i> strains, as long as they contain antigen H:h (in the antigen combination H:e,h) in slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin, 22.06.2016

 place, date



 signature Head of Quality Control