

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

Product	Anti-Salmonella H:h
Catalogue no.	TR1409
Package size	1 ml
4850619	7240920
Expiry date	2023-08
Intended use	Anti-Salmonella H:h records all <i>Salmonella</i> strains, as long as they contain antigen H:h (in the antigen combination H:e,h) in slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin, 29.09.2020
place, date

J. Leutz
signature Head of Quality Control