

## Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
 ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485


sifin diagnostics gmbh

Berliner Allee 317-321  
 13088 Berlin, Germany

T: +49 30 927 030-0  
 F: +49 30 927 030-30  
 info@sifin.de  
 www.sifin.de

<b>Product</b>	<b>Anti-Salmonella H:r</b>
Catalogue no.	TR 1416
Package size	1 ml
<b>Lot no.</b>	<b>149 09 16</b>
<b>Expiry date</b>	<b>2019-03</b>
Intended use	Anti-Salmonella H:r records all <i>Salmonella</i> strains, as long as they contain antigen H:r or H:r(i) in slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin, 24.10.2016  
 \_\_\_\_\_  
 place, date

  
 \_\_\_\_\_  
 signature Head of Quality Control