

## Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321  
13088 Berlin, Germany

T: +49 30 927 030-0  
F: +49 30 927 030-30  
info@sifin.de  
www.sifin.de

|                    |  |
|--------------------|--|
| <b>Product</b>     | <b>Anti-Salmonella H:r</b>   |
| Catalogue no.      | TR 1416  |
| Package size       | 1 ml   |
| <b>0461118</b>     | <b>0461118</b>   |
| <b>Expiry date</b> | <b>2021-05</b>   |
| Intended use       | Anti-Salmonella H:r records all <i>Salmonella</i> strains, as long as they contain antigen H:r or H:r(i) in slide agglutination. |
| Specificity        | No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens.                                |
| Preservatives      | Sodium azide 0.9 mg/ml   |
| Storage            | 2 °C...8 °C  |

Berlin,  
place, date

25.07.2019

  
signature Head of Quality Control