

sifin

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321
13088 Berlin, Germany

T: +49 30 927 030-0


F: +49 30 927 030-30

info@sifin.de

www.sifin.de

Product	Anti-Salmonella H:6
Catalogue no.	TR 1435
Package size	1 ml
Lot no.	199 10 14
Expiry date	2017-10
Intended use	Anti-Salmonella H:6 records all <i>Salmonella</i> strains, as long as they contain antigen H:1,6 in slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C... 8 °C

Berlin, 01.12.14
place, date


signature Head of Quality Control